Read instructions carefully before starting test

INTENDED USE

In an AOAC Research Institute Performance TestedSM Method study, ANSR for *Listeria* has been found to be an effective procedure for detection of *Listeria* spp. in deli turkey, hot dogs, fermented sausage, pasteurized milk, Mexican-style cheese, ice cream, smoked salmon, lettuce, cantaloupe, guacamole, pasteurized liquid egg and sponge or swab samples from stainless steel, plastic, ceramic, sealed concrete and rubber environmental surfaces.

Although the ANSR test chemistry detects *L. grayi*, detection of this species in enriched samples will be variable depending on the ability of individual strains to grow under the enrichment conditions used in the method.

ASSAY PRINCIPLES

ANSR for *Listeria* is an isothermal, amplified nucleic acid assay. The ANSR for *Listeria* method is based on nicking enzyme amplification reaction (NEARTM) technology preceded by the reverse transcription of 23S ribosomal RNA. Target complementary DNA is amplified through a mechanism of polymerization from the ends of nicks created in double-stranded DNA by the action of a specific endonuclease. Amplified target sequences are detected in real time using fluorescent molecular beacon probes.

A 2-stage lysis reaction is performed, first at $37 \pm 2^\circ\text{C}$ for 10 minutes, then at $80 \pm 2^\circ\text{C}$ for 20 minutes. Next, a portion of the lysed sample is transferred to a strip tube containing lyophilized ANSR reagents. The tubes are sealed and incubated at $56 \pm 2^\circ\text{C}$ on the ANSR reader. Results are generated by the reader and displayed in the ANSR software within 18 minutes. Positive results may be confirmed from the enrichment cultures following standard procedures. Each tube of ANSR reagents contains an internal positive control, ensuring that the reagents are functioning properly.

INTENDED USER

The ANSR for *Listeria* test is designed for use by personnel with appropriate training in microbiology. Training in the use of the ANSR test system is available through Neogen.

MATERIALS PROVIDED

1. 12 strips of 8 cluster tubes, 1.2 mL
2. 12 strips of 8 reaction tubes, 200 μ L, containing lyophilized ANSR for *Listeria* reagents in 2 sealed foil pouches with desiccant pack
3. 12 strips of 8 permanent caps for reaction tubes
4. 1 bottle lysis reagent suspension buffer, 60 mL
5. 3 vials containing lyophilized lysis reagents
6. 1 kit insert

EQUIPMENT REQUIRED

1. ANSR reader (Neogen item 9828)
2. Computer and software for connection to ANSR reader (Neogen item 9832)
3. 1 dual heater block with aluminum block inserts for 1.2 mL cluster tubes, $80 \pm 2^\circ\text{C}$ and $37 \pm 2^\circ\text{C}$ (Neogen items 9386-48D, 9829-48) **OR** two single heater blocks with aluminum block inserts for 1.2 mL cluster tubes, $80 \pm 2^\circ\text{C}$ and $37 \pm 2^\circ\text{C}$ (Neogen items 9386-48D, 9829-48)
4. Pipettor, 100–1000 μ L (Neogen item 9463)
5. Pipette tip rack, 100–1000 μ L, sterile, 96 tips (Neogen item 9487)
6. Pipettor, 10–100 μ L, 8-channel (Neogen item 9388)
7. Pipette tips, 100 μ L, sterile, filtered, 96 tips (Neogen item 9389)
8. Vortex mixer, adjustable speed (Neogen item 9494)
9. Stomacher (optional)
10. 3 thermometers, traceable (Neogen item 9518)
11. Timer, 3-channel (Neogen item 9426)
12. Optional-for use-heater block with 0.2 mL reaction tube aluminum block insert, $56 \pm 2^\circ\text{C}$ (Neogen items 9386-48D, 9829-48)
13. Webcam (Neogen item WEBCAM)
14. ANSR Ethernet cable (Neogen item 9835)
15. 10 mL pipette pump (Neogen item 9277)
16. Pipettes, sterile serological (Neogen item 8686)
17. 40-slot, 20 mm test tube rack, autoclavable (Neogen item 9553)
18. Pipettor, 20–200 mL (Neogen item 9276)

OTHER MATERIALS REQUIRED

1. Stomacher-type bags for sample enrichment. Filtered bags are recommended (Neogen item 6827)
2. Graduated cylinder, 250 mL (Neogen item 9368)

MEDIA ENRICHMENT BROTH REQUIRED

1. ANSR for *Listeria* Enrichment Broth 1 (Neogen item 9814) — for environmental samples
2. LESS Medium, 500 g bottle (Neogen item 9790A) — for food samples except meats
3. Reveal Buffered *Listeria* Enrichment Broth (BLEB), 500 g bottle (Neogen item 9786) — for meat samples. Also used for confirmation, see page 6.
4. LESS Plus Medium, 500 g bottle (Neogen item 9880) — for select matrices.
5. 1 L purified water
6. Phosphate Buffered Saline (if testing smoked salmon or guacamole) (Neogen item 8425)

STORAGE

Store ANSR reagents at 2–8°C. After removing reaction tubes from the foil pouch, promptly re-seal the pouch. Leave the desiccant pack in the pouch at all times.

PRECAUTIONS

1. Use good microbiology laboratory practices.
2. **Dispose of used pipette tips in a covered container containing a fresh solution of 10% bleach. The 10% bleach solution should be made fresh each day. Use 9 parts water with 1 part commercially available household-strength bleach to make the**

10% bleach solution. Undiluted stock solutions of bleach should be used within 30 days after opening.

3. **Discard bleach solution and tips as regular waste at the end of each day.**
4. *Listeria monocytogenes* is a known hazard to pregnant women and immunocompromised individuals. Consult with your facility safety director for specific instructions.
5. Do not use reagents beyond the expiration date.
6. Enrichment media should be prewarmed to incubation temperature (all protocols) before use.
7. Use of enrichment media and incubation times or temperatures other than those specified may lead to erroneous results.
8. Remove reaction tubes from the foil pouch just before use and keep covered until heating process begins. Reseal the pouch containing the remaining reaction tubes to **avoid prolonged exposure to light**. More than 15 minutes of total exposure time may lead to erroneous results.
9. Do **NOT**, under any circumstance, remove caps from reaction tubes **AFTER** the assay has been started. This is essential in order to prevent accidental contamination of the environment with amplification products.
10. Exercise care in all pipetting steps to avoid cross-contamination of samples.
11. Complete all assay steps in sequence, avoiding delays between steps.
12. Tap reaction tubes on bench top to make sure lyophilized reagents are at the bottom of the tube prior to adding lysed sample.

PREPARATION OF ENRICHMENT BROTH

Prepare all media according to manufacturer's instructions.

NOTE: Use rehydrated media the same day as prepared.

SAMPLE ENRICHMENT

Environmental samples

1. Place the sponge or swab sample in a Stomacher-type bag or tube.
2. Add an appropriate amount of **ANSR for *Listeria* Enrichment Broth 1** (Neogen item 9814) that has been prewarmed to $36 \pm 1^\circ\text{C}$. For sponge samples, an appropriate amount typically is 100–200 mL. For swab samples, it typically is 10 mL.
3. Incubate the culture at $36 \pm 1^\circ\text{C}$ for **16–24 hours**.

Food samples

1. Weigh out 25 g of sample in a Stomacher-type bag.
2. Add 225 mL LESS Medium prewarmed to $36 \pm 1^\circ\text{C}$ to the bag. Homogenize (Stomacher, etc.) the sample as appropriate for the sample type.
3. Incubate the culture at $36 \pm 1^\circ\text{C}$ for **16–24 hours**. For Mexican-style cheese, milk and smoked salmon, incubate **22–26 hours**.

NOTE: For larger sample sizes, maintain a 9:1 ratio of enrichment broth to sample.

Meat samples

1. Weigh out 125 g of sample in a Stomacher-type bag.
2. Add 1125 mL of Reveal BLEB prewarmed to $36 \pm 1^\circ\text{C}$ to the bag. Homogenize (Stomacher, etc.) the sample as appropriate for the sample type.
3. Incubate the culture at $36 \pm 1^\circ\text{C}$ for **22–26 hours**.

ALTERNATIVE SAMPLE ENRICHMENT – For select matrices

Environmental Samples

1. Place the sponge or swab sample in a Stomacher-type bag or tube.
2. Add an appropriate amount of LESS Plus Medium (Neogen item 9880) prewarmed to $36 \pm 1^\circ\text{C}$ to the bag or tube. 100 mL and 10 mL of medium are recommended for sponge and swab samples, respectively.
3. Incubate the culture at $36 \pm 1^\circ\text{C}$ for **16–24 hours**.

Cantaloupe

1. Weigh out 25 g of sample in a Stomacher-type bag.
2. Add 225 mL LESS Plus Medium (Neogen item 9880) prewarmed to $36 \pm 1^\circ\text{C}$ to the bag. Homogenize (Stomacher or equivalent) the sample as appropriate for the sample type.
3. Incubate the culture at $36 \pm 1^\circ\text{C}$ for **16–24 hours**.

Hot dogs and Mexican-style cheese

1. Weigh out 125 g of sample in a Stomacher-type bag.
2. Add 375 mL LESS Plus Medium (Neogen item 9880) prewarmed to $36 \pm 1^\circ\text{C}$ to the bag. Homogenize (Stomacher or equivalent) the sample as appropriate for the sample type.
3. Incubate the culture at $36 \pm 1^\circ\text{C}$ for **22–26 hours**.

LYSIS REAGENT SOLUTION PREPARATION

1. Reconstitute 1 vial of lyophilized lysis reagents with 18 mL of lysis reagent suspension buffer by adding the buffer to the reagent vial. Swirl gently to mix.
NOTE: 1 vial of lysis reagents is enough for approximately 32 samples. Prepared lysis reagent solution can be stored at $2\text{--}8^\circ\text{C}$ for up to 30 days.

ANSR TEST PROCEDURE

Prior to starting the assay:

1. Preheat one lysis heater block to $80 \pm 2^\circ\text{C}$. Preheat the second lysis heater block to $37 \pm 2^\circ\text{C}$. If using the optional single heater, preheat to $56 \pm 2^\circ\text{C}$. Use the thermometer for the temperature reading.
2. Remove the foil pouch containing the reaction tubes from the refrigerator and **allow them to warm at room temperature for 15 minutes**. To avoid excess light exposure, leave reaction tubes in foil pouch until they are needed.

NOTE: The lysis reagent solution should also be allowed to be at room temperature for 15 minutes before use.

3. Connect the ANSR reader to the computer via USB or Ethernet and turn the computer on.
4. Turn on the ANSR reader. The reader will preheat to $56 \pm 2^\circ\text{C}$.
5. Start the ANSR software and click the Connect button. Input sample IDs, lot number and user information.

NOTE: For instructions on using the reader and software, see the user guide that came with the ANSR reader.

ASSAY PROCEDURE

NOTE: Certain sample types require dilution of the enrichment culture prior to the assay. For smoked salmon, guacamole, pasteurized liquid egg and meat samples, dilute the enrichment culture 1:10 in Phosphate Buffered Saline (PBS) and use 50 μ L of the dilution for the ANSR assay.

1. Add 50 μ L enrichment culture to a 1.2 mL cluster tube(s) using a pipettor with 100 μ L filtered tips.

NOTE: Cluster tubes may be pulled apart to provide the number of tubes needed.

2. Add 450 μ L lysis reagent solution to each cluster tube(s) containing culture. Be sure to switch pipette tips between samples.

NOTE: Return the lysis reagent solution to the refrigerator after use (within 1 hour).

3. Incubate the cluster tube(s) at $37 \pm 2^\circ\text{C}$ for **10 minutes**.

4. Immediately transfer the cluster tube(s) to the $80 \pm 2^\circ\text{C}$ heater block and incubate for **20 minutes**.

NOTE: The $80 \pm 2^\circ\text{C}$ incubation time may be extended to a total of **60 minutes** for the purpose of managing staggered assay start times.

5. For **3 to 5 minutes** before the end of the lysis step, preheat the ANSR reagents to $56 \pm 2^\circ\text{C}$ by placing the reaction tubes in the ANSR reader. Optional: A separate heat block can be used. It should be heated to $56 \pm 2^\circ\text{C}$.

NOTE: The strip of reaction tubes may be cut to provide the number of tubes needed.

Keep all unused tubes in the sealed foil pack. Ensure the pellet in the reaction tube(s) is at the bottom by tapping the tubes gently on the bench top.

6. After the completion of the 20–60 minute lysis incubation, remove and discard the caps from the reaction tube(s) in the ANSR reader.

IMPORTANT: Proceed with steps 7–9 without delay. **The transfer of the sample from the lysis tubes to the reaction tubes should be completed within 1 minute.**

7. Using an 8-channel pipette and 100 μ L filter tips, carefully transfer 50 μ L from the top third of the lysed sample(s) in the cluster tube(s) to the reaction tube(s). Debris may accumulate at the bottom of the lysis tube(s) that will interfere with assay performance. Avoid transfer of debris by aspirating from the top third of the lysis tube(s). Do not prime the pipette tips and do not mix before aspirating.

Place the provided permanent cap(s) on the reaction tube(s).

NOTE: Lysed sample may be transferred from the same cluster tube a maximum of 3 times.

8. Remove the strips(s) of tubes from the reader (or $56 \pm 2^\circ\text{C}$ heat block if one was used) and vortex briefly (about 2 seconds) at medium/high speed (2500 RPM). Quickly visually check each reaction tube to assure that no bubbles are on the bottom or in the middle of the sample. A quick tap of the tubes should release any bubbles from the bottom or middle to the top. Then place into the reader without delay. Close the reader's lid.

NOTE: The reader will not provide accurate results if the lid is open. Keep the lid closed at all times while the assay is running. Contamination may occur if the permanent caps are not placed on the reaction tubes and/or if the permanent caps are removed.

9. Click **START** in the ANSR software to begin the **18 minute** assay.
10. Results will be displayed as positive, negative or invalid once the assay is finished. If the result displays an invalid, the test must be repeated. A single repeat test can be run starting at step 6. Alternatively, start over from step 1 if lysis from step 6 has been at $80 \pm 2^\circ\text{C}$ for more than 60 minutes.

INTERPRETATION OF RESULTS

The ANSR software will indicate the test results as positive or negative for the presence of *Listeria* spp. in the enriched sample. In addition, the real-time fluorescence curve generated from the assay can be viewed.

CONFIRMATION

Samples producing positive ANSR results may be confirmed by streaking the enrichment cultures to *Listeria* selective/differential agar media and continuing with identification of presumptive colonies following standard methods ^[1,2].

For certain sample types (e.g., processed meats, milk, cheese, ice cream, lettuce and guacamole), a secondary broth enrichment may be necessary prior to plating to reduce background flora. Transfer 0.1 mL of the primary enrichment culture to a tube containing either 10 mL of Buffered *Listeria* Enrichment BrothTM (BLEB) or, in the case of meat commodities, 10 mL of MOPS-BLEB. Incubate either the BLEB culture for **26 ± 2 hours** at 35 ± 2°C, or the MOPS-BLEB culture for **18-24 hours** at 35 ± 2°C, then proceed with streaking to selective/differential agar media and subsequent identification per standard methods.

DISPOSAL

Enrichment cultures and used lysis tubes should be disposed of as biohazard waste. The preferred method of treatment for biohazard waste is autoclaving. Items that **CANNOT** be autoclaved, such as used reaction tubes and caps, should be immersed in a fresh 10% bleach solution that is made daily. Consult with the safety advisor of your facility for detailed instructions.

Do **NOT** remove the permanent caps, for **ANY** reason, from the ANSR reaction tubes once the assay has started, even when disposing of them. Reaction tubes can be disposed of as non-biohazardous waste. It is recommended that they be placed in a sealable container or plastic bag with 10% bleach and immediately disposed of to protect against accidental opening.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

Neogen ANSR Molecular Diagnostic for Foodborne Pathogen Detection – limited use label license

SYTO® 82

SYTO® 82 contained within this product is provided under an intellectual property license from Life Technologies Corporation, Eugene, OR, and may be used for *in vitro* detection and analysis of (i) food, feeds and beverages, including nutraceuticals, (ii) ingredients for food, feeds and beverages, (iii) process samples from food, feed and beverage preparation, distribution and delivery, and (iv) water from any source for human consumption, all for the purpose of safety and quality assurance. The buyer must not sell or otherwise transfer this product or its components for any other use, including but not limited to: human *in vitro*, veterinary, identity or paternity testing, forensics, or *in vivo* detection of nucleic acid sequences in living beings, or cells. For information on purchasing a license for SYTO 82 for purposes other than food, beverage and water safety, and quality assurance, contact Life Technologies Corporation at outlicensing@lifetech.com.

Molecular beacon probes

One or more molecular beacon probes contained within this product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for tests on food products, feeds, beverages and water.

NEAR™ Technology

This product utilizes the patented NEAR isothermal technology and is sold under license from Ionian Technologies, San Diego, CA, and may be used under Ionian Technologies patent rights only for tests on food, beverage and water safety.

References

- [1] USDA-FSIS (2011) Microbiology Laboratory Guidebook, chapter 8.07
<http://www.fsis.usda.gov/wps/wcm/connect/1710bee8-76b9-4e6c-92fc-fdc290dbfa92/MLG-8.pdf?MOD=AJPERES>
- [2] USDA-FDA (2013) Bacteriological Analytical Manual, chapter 10
<http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/ucm071400.htm>

PRODUCTS AVAILABLE

- 9465 **Environmental sampling sponges** (with Dey Engley (DE) neutralizing broth) – pack of 100
- 36002 **Sampling sponges, pre-moistened, HiCap** – pack of 100
- 9822 **ANSR for *E. coli* O157:H7** – up to 96 tests
- 9824 **ANSR for *Listeria monocytogenes*** – up to 96 tests
- 9870 **ANSR for *Salmonella*** – up to 96 tests
- 9871 **ANSR for *Listeria*** – up to 96 tests
- 9872 **ANSR for *Campylobacter*** – up to 96 tests
- 9873 **ANSR for *Listeria* Right Now™** – up to 96 tests
- 9837 **Complete ANSR Equipment Kit** – Reader, computer, pipette set, dry bath heaters and inserts, thermometers, vortex, time, rack

**North America****Neogen Headquarters**

800/234-5333 (USA/Canada)
foodsafety@neogen.com
foodsafety.neogen.com

Europe, Middle East and Africa**Neogen Europe**

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China**Neogen Bio-Scientific Technology**

+86 21 6271 7013
info@neogenchina.com.cn
www.neogenchina.com.cn

India**Neogen Food and Animal Security**

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