



soleris

E. coli Vial

Product Number: S2-EC



Pictured: Soleris® vial uninoculated (left) and inoculated vial (right).

Introduction

The *E. coli* vial (S2-EC) is used to detect *Escherichia coli* species. The S2-EC vial contains a broth enrichment medium with a fluorescent indicator. The vial has an assay time of 18 hours for most applications. The reading frame contains a soft agar plug which reduces product interference with the detection of microbial growth. As *E. coli* organism grows in the broth medium, the utilization of MUG (4-methylumbelliferyl-beta-D-glucuronide) by *E. coli* enzymatic reaction produced diffuses into the soft agar plug. The fluorescent end product called methylumbelliferone is read by the optical sensors in the instrument.

Materials Required

1. S2-EC, *E. coli* (S2-EC) vial
2. Tryptic Soy Broth, (BLX-TSB90)
 - a. If required, use a designated neutralization broth, such as D/E Neutralizer, TAT Broth, Modified Letheen Broth, etc.
3. Soleris® *E. coli* MUG supplement (S2-ECM)
4. Kovac's Indole reagent (KO-122)
5. Test tube for confirmation

Dependent on Sample Tested

1. Sterile 1 N to 5 N sodium hydroxide (NaOH) and/or hydrochloric acid (HCl)
2. pH meter or pH paper

Vial Specifications

1. Vial pH is 7.4 ± 0.2
2. Vial sample capacity: 0.1 mL

Sample Preparation

1. For USP testing, perform 1:10 dilution by adding 10 g of sample in 90 mL of Tryptic Soy Broth or designated neutralization broth.
 - a. Check pH and adjust, if necessary, to 7.0 ± 1.0 .
2. Incubate for 18–24 hours at 35°C.

Inoculation of Vial

1. Transfer 0.12 mL of the Soleris *E. coli* MUG supplement (S2-ECM) to the S2-EC vial.
2. Cap the vial and gently invert 3 times to mix the sample.
3. Transfer 0.1 mL of the incubated enrichment to the S2-EC vial.
4. Cap the vial and gently invert 3 times to mix the sample. Keep the cap tight.



5. Insert the vial into the Soleris® Next Generation instrument set at 35°C and run for the preprogrammed test duration. It is not recommended to adjust the parameters without consulting NEOGEN® Technical Services.
6. If detection occurs, perform the Kovac's Indole confirmation test.

Algorithm Utilized

Test	Threshold	Skip	Shuteye	Test Duration	Temperature
S2-EC	8	2	20	24 Hours	35°C

E. coli Confirmation Step (Kovac's Indole)

1. Remove the S2-EC vial positive (detecting) vial from the instrument.
2. Transfer 0.5–1.0 mL to a test tube.
3. Add a few drops (5–8) of Kovac's reagent (KO-122)
 - a. Do not add Kovac's directly to the S2-EC vial, as it will destroy all viable organisms for identification.
4. Appearance of a bright pink ring at the meniscus of the broth indicates the presumptive presence of *E. coli*.
 - a. If negative, the reagent layer will retain the original yellow color.
5. Presumptive positive samples should be sent out for identification.

Alternative Testing Using MacConkey Broth

1. Follow appropriate neutralization scheme as recommended in sample neutralization procedures. Adjust the pH to 7.0 ± 1.0 , if necessary.
2. Incubate the sample at $35 \pm 0.5^\circ\text{C}$ for 18–24 hours.
 - a. If testing for *E. coli*, transfer 0.1 mL of overnight enriched sample to 10 mL of MacConkey Broth and incubate for 18–24 hours at 35°C.
3. Supplement the appropriate vial:
 - a. S2-EC — add 0.120 mL of S2-ECM Supplement to the vial.
4. Transfer 0.100 mL of the homogenized, overnight enriched sample (for S2-EC transfer from enriched MacConkey Broth) into the vial.
 - a. Mix the vial contents and insert them into the instrument at the recommended testing parameters.
5. If positive results occur:
 - a. In the vial, perform the confirmation method described in the vial technical information sheet.

Disclaimers

Information provided is based on validation procedures that NEOGEN performed in NEOGEN Laboratories. Deviation from procedures is possible, but should be discussed with NEOGEN Technical Services.

Samples may need to be pH adjusted for all vials.

Appearance of the vials should be inspected prior to use.

If shuteye detections are observed, the threshold may need to be adjusted based on the product matrix. Certain product matrices may require parameter adjustments, including increased test duration. For more information contact NEOGEN Technical Services.



Pictured: Confirmed presumptive positive (left) and confirmed negative (right).