

Ag-Tek® Veterinary Apparel

Designed for veterinarian use when performing cesarean sections and other large animal procedures, Ag-Tek veterinary apparel is made from a urethane coated nylon oxford material. The garments are navy blue in color, highly water resistant, breathable and machine washable. Except for the Shoulder Protector, BodyGuard and Vet Apron, all garments are available in four sizes: S, M, L and XL. Custom sizes of all garments may be special ordered. Please specify size when ordering, by adding S, M, L or XL after the item number (for example, VA100-S).

Apparel Measurements

Tops

Size	Chest	Arms
SMALL	46"	18"
MEDIUM	49"	18"
LARGE	54"	18"
X-LARGE	55"	20"
XX-LARGE	64"	25"

Bottoms

Size	Waist	Inseam
SMALL	37"	30"
MEDIUM	42"	31"
LARGE	44"	32"
X-LARGE	48"	34"
XX-LARGE	49"	34"

Cleaning Instructions

- Shake out after use
- Remove any excess dirt and grime
- Spray with hose and pre-soak
- Wash - normal cycle in COLD water
- Shake out - DO NOT DRY IN DRYER - line dry only

Vet Jacket

- Jacket features attached hood

Item No.	Description
VA300	Vet Jacket (Please specify size when ordering: S, M, L or XL)

OB Jacket

Item No.	Description
VA301	OB Jacket (Please specify size when ordering: S, M, L or XL)

OB Suit

- Pants have an elastic waistband

Item No.	Description
VA101	Top only - specify size: S, M, L or XL
VA102	Bottoms only - specify size: S, M, L or XL

Vet Apron

- Extra long ties

Item No.	Description
VA400	Vet Apron - One size

OB Gown

- Velcro fasteners in back

Item No.	Description
VA500	OB Gown (Please specify size when ordering: S, M, L or XL)

BodyGuard™

- Buckle with adjustable strap
- May be worn on either side

Item No.	Description
VA800	BodyGuard (Please specify size when ordering: S or L)



OB Coverall

- Two-way zipper

Size	Chest	Waist	Inseam
SMALL	42"	42"	20"
MEDIUM	44"	44"	23"
LARGE	48"	48"	26"
X-LARGE	50"	50"	29"

Item No.	Description
VA700	OB Coverall (Please specify size when ordering: S, M, L or XL)

Vet Overall

- Velcro fasteners on straps

Size	Waist	Inseam
SMALL	36"	31"
MEDIUM	39"	32"
LARGE	41"	33"
X-LARGE	47"	33½"

Item No.	Description
VA600	Vet Overall (Please specify size when ordering: S, M, L or XL)

Shoulder Protector

- Buckle with adjustable strap
- May be worn on either arm

Item No.	Description
VA900	Shoulder Protector (Please specify size when ordering: S or L)



OB Coverall



Vet Overall



Shoulder Protector



PolyGown™



PolyApron™



Veterinary Gown

Disposable Apparel

PolyGown™

Manufactured of 1.3 mil lightweight embossed polyethylene film, the Ag-Tek PolyGown is designed for veterinary and livestock applications where minimal protection is needed. PolyGown is packaged in a convenient dispenser box for easy removal and is available in three sizes: M, L and XL.

Item No.	Description
DG100-M	PolyGown, size M - box of 20
DG100-L	PolyGown, size L - box of 15
DG100-XL	PolyGown, size XL - box of 10

PolyApron™

Designed for durable protection of clothing when treating or grooming large and small animals, the disposable Ag-Tek PolyApron is made from blue heavyweight 4.5 mil polyethylene film. Available in one size (28" x 46"), the apron has extra-long ties for ease of tying in front. PolyApron is packaged and sealed in a polyethylene sleeve containing 10 aprons.

Item No.	Description
A450	PolyApron, one size - pack of 10

Veterinary Gown

The Ag-Tek Veterinary Gown is manufactured of durable 1.2 mil polyethylene film and is designed for those veterinary procedures where disposable full-body protection is needed. The Veterinary Gown is 61" in length, has a slit neck for ease in donning and comfort pleats in the back to facilitate ease of movement. The Veterinary Gown is packaged in a heat-sealed bag containing 25 gowns. One size.

Item No.	Description
VG100	Veterinary Gown, one size - pack of 25