



# Ideal<sup>®</sup> Disposables

The *Ideal* Choice for Disposables



944 Nandino Blvd. • Lexington, KY 40511  
859/254-1221 or 800/621-8829 • fax: 800/621-8829  
inform@neogen.com • www.neogen.com




Trusted. Valued. Supported.  
ANIMAL SAFETY SOLUTIONS

# Poly Hub Disposable Needles

# Disposables Retail Packs

## Poly Hub Features - Cartridge Packed:

- Clear hubs for easy viewing of flashback
- Metal insert increases needle strength for animal health applications
- Tri-beveled, anti-coring tips are computer-ground for consistent sharpness
- Anti-friction coating for minimal penetration resistance and animal discomfort
- Stainless steel cannula
- Sterile, nonpyrogenic, autoclavable components
- Tamper evident cartridges are safe, convenient, easily disposed of and transportable
- Lot number is printed on each cartridge for easy tracking
- ISO color-coded cartridges for easy gauge identification



Needles are available in boxes of 100; master cases contain 10 boxes.

Ask about additional size and length availability.

Item No.	Needle Gauge/Length	Item No.	Needle Gauge/Length
9319	16 G x 1/2"	9365	20 G x 3/4"
9325	16 G x 1"	9370	20 G x 1"
9330	16 G x 1 1/2"	9375	20 G x 1 1/2"
9340	18 G x 3/4"	9390	22 G x 3/4"
9345	18 G x 1"	9395	22 G x 1"
9350	18 G x 1 1/2"	9399	25 G x 5/8"
9360	20 G x 1/2"		

## Poly Hub Features - Soft Packed:

- Ultra-sharp needle
- Clear polypropylene hubs for easy viewing of flashback
- Paper-to-plastic packaging is tamper evident and ensures sterility

Needles are available in boxes of 100; master cases contain 10 boxes.

Ask about additional size and length availability.

Item No.	Needle Gauge/Length	Item No.	Needle Gauge/Length
9345K	18 G x 1"	9390K	22 G x 3/4"
9350K	18 G x 1 1/2"	9395K	22 G x 1"
9370K	20 G x 1"	9393K	22 G x 1 1/2"
9375K	20 G x 1 1/2"	9398K	23 G x 1"
9380K	21 G x 1"	9399K	25 G x 5/8"



## Aluminum Hub Disposable Needles

Item No.	Needle Gauge/Length	No. Per Package
9323	16 G x 1"	5
9329	16 G x 1"	25
9328	16 G x 1 1/2"	5
9342	18 G x 5/8"	5
9347	18 G x 1"	5
9348	18 G x 1"	25

## Poly Hub Disposable Needles

Item No.	Needle Gauge/Length	No. Per Package
9318	16 G x 1/2"	5
9326	16 G x 1"	5
9327	16 G x 1"	25
9331	16 G x 1 1/2"	5
9341	18 G x 3/4"	5
9346	18 G x 1"	5
9355	18 G x 1"	25
9351	18 G x 1 1/2"	5
9352	18 G x 1 1/2"	25
9361	20 G x 1/2"	5
9362	20 G x 1/2"	25
9367	20 G x 3/4"	5
9366	20 G x 3/4"	25
9372	20 G x 1"	5
9373	20 G x 1 1/2"	5
9391	22 G x 3/4"	5
9392	22 G x 3/4"	25
9397	22 G x 1"	5
9396	22 G x 1"	25

## Disposable Syringes, Luer Slip

Item No.	Size	No. Per Package
9276	3 cc	6
9271	6 cc	6
9272	12 cc	4
9273	20 cc	4
9274	35 cc	2
9275	60 cc	2

## Disposable Syringes, Luer Lock

Item No.	Size	No. Per Package
9263	3 cc	6
9264	6 cc	6
9265	12 cc	4
9266	20 cc	4
9267	35 cc	2
9268	60 cc	2

## Disposable Syringes, Luer Lock with Poly Hub Needles

Item No.	Needle Gauge/Length	No. Per Package
9184	3 cc, 22 G x 3/4"	3
9185	3 cc, 22 G x 1"	3
9187	12 cc, 18 G x 1"	3
9188	12 cc, 20 G x 1"	3



Ideal utilizes Plan-O-Gram retail display software.

## Ask about our other needle & syringe products:

- TU-MasterLine™ Products
- MegaShot™ 50cc Syringe
- Pro-Shot™ & Pro-Shot™ II Syringe
- Stainless Steel Durable Needles

### ISO Needle Gauge Color Code Chart

14 G	16 G	18 G	20 G	21 G	22 G	23 G	25 G
Green	White	Pink	Yellow	Dark Green	Grey	Blue	Orange

# D3 Detectable Needles™

# Aluminum Hub Disposable Needles

## Durable. Detectable. Dependable

### D3 Detectable Features: Detectability

Ideal has developed and added a detectable alloy to the cannula of the D3 Needles. This patented alloy allows D3 Needles to be as much as 100% detectable in metal detection systems, while maintaining their stainless steel classification.

### Strength

D3 Needles have thick-walled cannulas, which are crimped and pressed into sturdy, plated brass hubs. These stronger sidewalls increase load strength. Bending and breaking is reduced more than 6-fold over conventional needles, and needle tips stay sharp longer. The D3 Needle is strong enough to stand up to current vaccination practices without breakage.


### Sharpness

D3 Needles are ultrasharp with three-point precision ground tips for sharper cuts and easier penetration. Sharpness increases the possible number of injections per needle, and reduces the risk of injection site abscesses.

### Identification

Each D3 Needle has its gauge size followed by the letter D stamped into the brass hub. This allows users to easily distinguish detectable needles from conventional needles.

D3 Needles are specifically designed for production animal use and processing plant detectability. In addition, they address packer HACCP and quality assurance programs.



Needles are available in boxes of 100; master cases contain 10 boxes.

Item No.	Needle Gauge/Length
9418AL	14 G x 1"
9419AL	14 G x 1½"
9415	16 D x ⅝"
9416	16 D x ¾"
9420	16 D x 1"
9417	16 D x 1½"

Item No.	Needle Gauge/Length
9422	18 D x ⅝"
9421	18 D x ¾"
9423	18 D x 1"
9424	18 D x 1½"
9412AL	20 G x ½"
9411AL	20 G x 1"



### Aluminum Hub Features:


- Aluminum hub for large animal applications requiring a durable needle
- Tri-beveled, anti-coring tips are computer-ground for consistent sharpness
- Anti-friction coating for minimal penetration resistance and animal discomfort
- Stainless steel cannula
- Sterile, nonpyrogenic, autoclavable components
- Tamper evident cartridges are safe, convenient, easily disposed of and transportable
- Lot number is printed on each cartridge for easy tracking
- ISO color-coded cartridge for easy gauge identification

Needles are available in boxes of 100; master cases contain 10 boxes.

Ask about additional size and length availability.

Item No.	Needle Gauge/Length
9400AL	14 G x 1" <span style="color: green;">■</span>
9401AL	14 G x 1½" <span style="color: green;">■</span>
9402AL	14 G x 2" <span style="color: green;">■</span>
9403AL	16 G x ⅝" <span style="color: white;">■</span>
9404AL	16 G x ¾" <span style="color: white;">■</span>
9405AL	16 G x 1" <span style="color: white;">■</span>
9406AL	16 G x 1½" <span style="color: white;">■</span>

Item No.	Needle Gauge/Length
9407AL	18 G x ⅝" <span style="color: pink;">■</span>
9408AL	18 G x ¾" <span style="color: pink;">■</span>
9409AL	18 G x 1" <span style="color: pink;">■</span>
9410AL	18 G x 1½" <span style="color: pink;">■</span>
9413AL	20 G x 1" <span style="color: yellow;">■</span>
9414AL	20 G x 1½" <span style="color: yellow;">■</span>



## Satisfaction Guaranteed

We stand behind the quality of our disposables. If for any reason you are not 100% satisfied, Ideal will provide a refund or exchange the product free of charge.

## How do we ensure that our quality guidelines are met?

Critical quality specifications, which are derived from established industry standards. In total, there are 15 quality control specifications for our needles. If any of the measurements fall outside these specifications, the lot is rejected, ensuring that our needles are of the highest quality.

# Premium Disposable Syringes & Combos

# Standard Disposable Syringes & Combos

## Premium Syringe Features:

- Clear polypropylene barrel for easy viewing of fluid levels
- Large, easy to read graduations
- Positive plunger stop allows for maximum filling & eliminates costly spills
- Medical grade silicon lubricate inside barrel for smooth, accurate delivery
- Serrated thumb press for positive grip

## Premium Needle Features:

- Ultra-sharp, tri-beveled, anti-coring tips
- Metal inserts to increase strength for animal health applications
- Anti-friction coating for minimal penetration resistance and animal discomfort
- Clear polypropylene hubs for easy viewing of flashback

## Premium Packaging Features:

- Hard Pack design for veterinary office or field use
- ISO color-coded cartridge caps for easy gauge identification
- Sonic welded end cap provides tamper evident seal and insures sterility
- Product information sticker on box and case box for easy inventory control

### 1cc Tuberculin Syringes (Box of 100; case of 10 boxes)

Item No.	Needle Gauge/Length	Style
8859	25 G x 5/8"	Luer slip

### 3cc Syringes (Box of 100; case of 10 boxes)

Item No.	Needle Gauge/Length	Style
8861	N/A	Luer slip
8875	N/A	Luer lock
8877	20 G x 1"	Luer lock
8863	20 G x 1"	Luer slip
8878	20 G x 1 1/2"	Luer lock
8879	20 G x 1 1/2"	Luer slip
8864	22 G x 3/4"	Luer slip
8865	22 G x 3/4"	Luer lock
8866	22 G x 1"	Luer slip
8868	22 G x 1"	Luer lock
8876	25 G x 5/8"	Luer slip
8874	25 G x 5/8"	Luer lock



### 6cc Syringes (Box of 50; case of 10 boxes)

Item No.	Needle Gauge/Length	Style
8881	N/A	Luer slip
8882	N/A	Luer lock

### 12cc Syringes (Box of 80; case of 6 boxes)

Item No.	Needle Gauge/Length	Style
8883	N/A	Luer slip
8885	N/A	Luer lock
8884	18 G x 1"	Luer slip
8886	18 G x 1"	Luer lock

### 20cc Syringes (Box of 50; case of 6 boxes)

Item No.	Needle Gauge/Length	Style
8889	N/A	Luer slip
8890	N/A	Luer lock

### 35cc Syringes (Box of 30; case of 6 boxes)

Item No.	Needle Gauge/Length	Style
8891	N/A	Luer slip
8892	N/A	Luer lock

### 60cc Syringes (Box of 20; case of 6 boxes)

Item No.	Needle Gauge/Length	Style
8893	N/A	Luer slip
8894	N/A	Luer lock
8895	N/A	Catheter tip

## Standard Syringe Features:

- Clear polypropylene barrel for easy viewing of fluid levels
- Large, easy to read graduations
- Positive plunger stop allows for maximum filling and eliminates costly spills
- Medical grade silicon lubricate inside barrel for smooth, accurate delivery
- Serrated thumb press for positive grip

## Standard Needle Features:

- Ultra-sharp, tri-beveled, anti-coring tips
- Translucent ISO color coded hubs for easy viewing of flashback and quick gauge identification in or out of the packaging
- Anti-friction coating for minimal penetration resistance and animal discomfort

## Standard Packaging Features:

- Paper-to-plastic packaging is tamper evident and ensures sterility
- Product information sticker on box and case box for easy inventory control
- Each syringe/combo has UPC, lot number, manufacture and expiration dates printed on package

### U-40 Insulin Syringes (Box of 100; case of 10 boxes)

Item No.	Needle Gauge/Length
9101	1 cc w/ 29 G x 1/2"
9102	1/2 cc w/ 29 G x 1/2"

### 1cc Tuberculin Syringes (Box of 100; case of 24 boxes)

Item No.	Needle Gauge/Length	Style
9279	25 G x 5/8"	Luer slip

### 3cc Syringes (Box of 100; case of 24 boxes)

Item No.	Needle Gauge/Length	Style
9100	N/A	Luer slip
9170	N/A	Luer lock
9280	18 G x 1 1/2"	Luer lock
9281	20 G x 1"	Luer lock
9200	20 G x 1"	Luer slip
9282	20 G x 1 1/2"	Luer lock
9210	22 G x 3/4"	Luer slip
9284	22 G x 3/4"	Luer lock
9215	22 G x 1"	Luer slip
9285	22 G x 1"	Luer lock
9216	22 G x 1 1/2"	Luer slip



### 6cc Syringes (Box of 100; case of 24 boxes)

Item No.	Needle Gauge/Length	Style
9110	N/A	Luer slip
9171	N/A	Luer lock

### 12cc Syringes (Box of 100; case of 12 boxes)

Item No.	Needle Gauge/Length	Style
9120	N/A	Luer slip
9172	N/A	Luer lock
9250	18 G x 1"	Luer slip
9287	18 G x 1"	Luer lock
9288	20 G x 1"	Luer lock
9289	20 G x 1 1/2"	Luer lock

### 20cc Syringes (Box of 50; case of 18 boxes)

Item No.	Needle Gauge/Length	Style
9130	N/A	Luer slip
9173	N/A	Luer lock

### 35cc Syringes (Box of 50; case of 6 boxes)

Item No.	Needle Gauge/Length	Style
9140	N/A	Luer slip
9174	N/A	Luer lock

### 60cc Syringes (Box of 25; case of 8 boxes)

Item No.	Needle Gauge/Length	Style
9150	N/A	Luer slip
9175	N/A	Luer lock