

## Ag-Tek® MaxiSleeve®

Disposable MaxiSleeve OB sleeves are made from a blended polymer film for an optimum strength to sensitivity ratio. Packaged in resealable, pressure-sealed bags (PolyPack).

### M100 (Yellow)

- 1.10 mil thickness
- 39" overall length
- Designed for extra sensitivity
- Also available sterile (SM100)

### M125 (Brown)

- 1.25 mil thickness
- 39" overall length
- Durable copolymer film
- Also available sterile (SM125) and in fingerless design (M125-F)



Item No.	Sleeve Description
M100	MaxiSleeve - Yellow, 1.10 mil - pack of 100
M125	MaxiSleeve - Brown, 1.25 mil - pack of 100

## Ag-Tek® PolySleeve®

Disposable, polyethylene Poly-Sleeve OB sleeves are packaged in resealable, pressure-sealed bags (PolyPack).

### P100 (Buff; Light Duty)

- 1.10 mil thickness
- 36" overall length
- Ideal for AI procedures
- Also available sterile (SP100)

### P125 (Blue; Heavy Duty)

- 1.25 mil thickness
- 36" overall length
- Excellent strength/sensitivity ratio
- Also available sterile (SP125) and fingerless design (P125-F)

### P125-S (Blue; Heavy Duty)

- 1.25 mil thickness
- 36" overall length
- Same quality as P125 with smaller hand design
- Also available sterile small hand (SP125-S)

Item No.	Sleeve Description
P100	PolySleeve - Buff, 1.10 mil - pack of 100
P125	PolySleeve - Blue, 1.25 mil - pack of 100
P125-S	PolySleeve Small Hand - Blue, 1.25 mil - pack of 100



## Ag-Tek® MaxiSleeve® Heavy Duty

The Heavy Duty OB sleeve is made for the heavy duty jobs, being much thicker than other sleeves offered on the market. This glove is the same shape and length as the original, brown MaxiSleeve. Packaged in a convenient dispenser box.

- Same shape as MaxiSleeve
- Same length - 39" - as MaxiSleeve
- Thicker than MaxiSleeve at 1.4 mil thickness



Item No.	Sleeve Description
M125-HD	MaxiSleeve Heavy Duty OB Sleeve - Black, 1.4 mil - pack of 100


## Ag-Tek® PolyPetite™

*The smallest hand of any OB sleeve on the market*

The PolyPetite OB sleeve is fitted more snugly through the hand, wrist and arm. Being the smallest OB sleeve offered on the market, it is ideal for users with small hands. It is 36" in length and pink in color. Packaged in a convenient dispenser box.

- 1.25 mil thickness
- 36" length



 A portion of sale proceeds will go to breast cancer research at the Markey Cancer Foundation.

Item No.	Sleeve Description
P125-P	PolyPetite OB Sleeve - Pink - box of 100

## Ag-Tek® Sterile OB Gloves

For use in large animal surgery where sterility is required. Sterile versions of our MaxiSleeve, PolyPetite, PolySleeve, and EquiSleeve shoulder length sleeves are available individually wrapped in packs of 12.

Item No.	Sterile Sleeve Description
S-307620	Genia Extra Protection - Clear, 1.0 mil - pack of 12
SM100	MaxiSleeve - Yellow, 1.10 mil - pack of 12
SM125	MaxiSleeve - Brown, 1.25 mil - pack of 12
SM125-F	MaxiSleeve Fingerless - Brown, 1.25 mil - pack of 12
SP100	Poly-Sleeve - Buff, 1.10 mil - pack of 12
SP125	Poly-Sleeve - Blue, 1.25 mil - pack of 12
SP125-F	Poly-Sleeve Fingerless - Blue, 1.25 mil - pack of 12
SP125-P	PolyPetite - Pink, 1.25 mil - pack of 12
SP125-S	Poly-Sleeve Small Hand - Blue, 1.25 mil - pack of 12
SES100	EquiSleeve - Clear, 1.0 mil - pack of 12



Individual samples available upon request

## Ag-Tek® BreederSleeve®

Affordable, quality shoulder length sleeves. Packaged in a convenient dispenser box.

- 1.0 mil thickness
- 36" overall length
- Ideally suited for AI and breeding applications



Item No.	Sleeve Description
BS100	BreederSleeve - Red, 1.0 mil - box of 100

## Ag-Tek® BreederSleeve® Plus

Economical extra long sleeves, packaged in a convenient dispenser box.

- 1.0 mil thickness
- 39" overall length
- Ideally suited for AI and breeding applications



Item No.	Sleeve Description
BSP100	BreederSleeve - Orange, 1.0 mil - box of 100

## Poly-Shield™

Economical alternative to the BreederSleeve, designed for use where minimal protection is required. Packaged in a convenient dispenser box.

- 0.8 mil thickness
- 34" overall length



Item No.	Sleeve Description
PS151	Poly-Shield - Purple, 0.8 mil - box of 200

## Ag-Tek® Fingerless OB Sleeves for maximum finger sensitivity

Fingerless sleeves offer a veterinarian or technician full shoulder length protection, while providing the maximum degree of finger sensitivity. Our fingerless 1.25 mil MaxiSleeve and PolySleeve gloves are designed for use bare fingered or with a latex or nitrile examination glove. Sleeves are packaged in a resealable, pressure-sealed bag (PolyPack).

Item No.	Sleeve Description
M125-F	MaxiSleeve Fingerless - Brown, 1.25 mil - pack of 100
P125-F	PolySleeve Fingerless - Blue, 1.25 mil - pack of 100



## Ag-Tek® EquiSleeve®

Disposable EquiSleeve OB sleeves are made of a thin clear 1.0 mil copolymer film, providing the ultimate in sensitivity and softness. Due to the unique copolymer film, the EquiSleeve has extremely smooth seams reducing irritation of the sensitive mucosal lining, while providing the strength necessary for the most difficult exams. Packaged in resealable, pressure sealed bags (PolyPack).

- Ultimate in sensitivity and softness
- 1.0 mil thickness
- 36" overall length
- Also available sterile (SES100)



Item No.	Sleeve Description
ES100	EquiSleeve - Clear, 1.0 mil - pack of 100

## PolySem® Polytactyl®

The ultimate in strength and touch sensitivity. Packaged in a convenient dispenser box.

- 1.30 mil thickness
- 36" overall length
- Provides excellent touch sensitivity and optimum protection



Item No.	Sleeve Description
VGT3090	Polytactyl - Yellow, 1.30 mil - box of 100

## PolySem® Special Deliverance®

Disposable Special Deliverance sleeves provide excellent strength and are ideally suited for most animal health applications. Packaged in a convenient dispenser box.

- 1.20 mil thickness
- 36" overall length
- Suitable for a variety of applications



Item No.	Sleeve Description
VGS3090	Special Deliverance - Orange, 1.2 mil - box of 100

## PolySem® Polyver®

Disposable Polysem Polyver sleeves offer the thickest protection available and optimum resistance to tearing. Packaged in a convenient dispenser box.

- 1.40 mil thickness
- 36" overall length
- Thickest glove available



Item No.	Sleeve Description
VGP3090	Polyver - Green, 1.4 mil - box of 80

## Genia® Extra Protection OB Sleeve

This OB sleeve provides extra protection by covering shoulder, chest and under-arm area. The sleeve is held firmly in place by use of support strap. Sleeve measures 48" from fingertip to support strap. Packaged in a convenient dispenser box.



Item No.	Sleeve Description
307620	Extra Protection OB Sleeve - Clear, 1.0 mil - box of 50



Individual samples available upon request

Neogen Corporation • 944 Nandino Blvd. • Lexington, KY 40511 • 800/621-8829 • 859/254-1221 • inform@neogen.com • www.neogen.com