



Neogen microplate reader (6704)

- Fully automatic
- Bichromatic optics
- Includes four filters: 405, 450, 492, 650 nm
- On board data reduction calculations include regressions, curve fitting and cutoff absorbance calculations with step-by-step user prompting
- Large, non-volatile memory stores at least 50 user tests
- Includes parallel and serial outputs
- NRTL and CE certified

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com