

Reveal 3-D

for Food Allergens



Reveal 3-D Technology

The UNIQUE hand-held device with 3-line readout: a control line confirms the method has been performed successfully and two further lines differentiate low and high levels of detection.

Simple tools for food allergen control

The Reveal 3-D Tests make it simple to screen for the presence of low levels of allergens in CIP rinses and environmental swabs virtually anywhere. Sample preparation and testing take less than 10 minutes. Each kit contains everything needed (including swabs) to perform 10 single use tests (100 test kits also available).

Ordering Information

Prod.#	Product description
902086G	Reveal 3-D Almond – 10 samples
902075M	Reveal 3-D Casein – 10 samples
902082Q	Reveal 3-D Egg – 10 samples
901031P	Reveal 3-D Gluten – 10 samples
902087E	Reveal 3-D Hazelnut – 10 samples
901041L	Reveal 3-D Peanut – 10 samples
902081S	Reveal 3-D Crustacea – 10 samples
902093K	Reveal 3-D Soy – 10 samples

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com