

Reveal®

For ASP (Domoic Acid)




Intended Use

Reveal for ASP (Amnesic Shellfish Poisoning) is a 10 minute lateral flow assay that detects domoic acid. Domoic acid is a marine biotoxin that is secreted by algal blooms and is found in shellfish. Once the toxin is present in the shellfish tissue it cannot be destroyed by heat during cooking. Reveal for ASP detect domoic acid at the (established) U.S. and E.U. maximum threshold of 20 ppm.

Product Specifications

Sensitivity:	20 ppm
Testing time:	10 minutes
Antibody cross-reactivity:	Domoic Acid
Storage:	Room temperature, 18–30°C
Tests per kit:	24

Results



Ordering Information

Prod.#	Product description
9560	Reveal for ASP
9563	Reveal for Shellfish Toxin start-up kit

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com