

Trusted Results From Your Trusted Partner



Next Generation Reader

for
Neogen's Lateral Flow Technology

AccuScan III removes subjectivity from reading Reveal® lateral flow devices, provides consistent interpretation, and the accurate records to back it up. Tests include pathogens, mycotoxins, allergens, dairy antibiotics, GMOs and ruminant by-products.



We stand behind our results.


Consistent, accurate, traceable results that are easily analyzed.

AccuScan III offers:


- » **Consistency.**
AccuScan III's digital reading eliminates subjectivity from sample to sample, and technician to technician.
- » **Permanence.**
By recording the sample ID, digital image of the results, time and date, AccuScan III eliminates manual recording, and provides easy reporting and permanent, traceable results.
- » **Simple PC interface.**
AccuScan III's test results can be exported as digital pictures, and to Microsoft Excel for powerful and permanent reporting, and tracking and trend analysis.
- » **Intuitive operation.**
AccuScan III's personal data assistant platform minimizes training needs, and its small and lightweight design allows for easy transportation between locations.
- » **Flexibility.**
AccuScan III is designed for use with all of Neogen's lateral flow products, including tests for pathogens, mycotoxins, food allergens, dairy antibiotics, GMOs, and ruminant by-products.
- » **Dependability.**
AccuScan III is from Neogen, a proven leader in food safety innovation that has been in business since 1982. We take pride in our unequaled, 24-hour customer service.



Optional Printer



Results screen



DataManager



800/234-5333 or 517/372-9200
620 Leshar Place, Lansing, MI 48912
foodsafety@neogen.com • www.neogen.com