

Reveal®

For Ruminant in MBM




Intended Use

Reveal for Ruminant in MBM is an immunochromatographic lateral flow assay used for the qualitative analysis of ruminant MBM in nonruminant animal meal. The test can detect as little as 2% ruminant MBM in nonruminant animal meal.

Product Specifications

Sensitivity:	2% in meat and bone meal
Testing time:	Extraction boiling step: 10 minutes; Test strip development: 10 minutes
Storage:	Room temperature, 18–30°C
Tests per kit:	25
Approvals:	AOAC-RI #070501

Results



Ordering Information

Prod.#	Product description
8105	Reveal for Ruminant in MBM - 25 samples
9590	Reveal AccuScan III Reader

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com