

Reveal®

For *Salmonella Enteritidis*




Intended Use

The Reveal for *Salmonella Enteritidis* (SE) Test System provides for the rapid recovery of *Salmonella Enteritidis* in food, layer or broiler house environments and animal feed, allowing detection and presumptive identification of the target organism.

Product Specifications

Sensitivity:	1–10 CFU/25 g
Testing time:	15 minutes (after 48 hour enrichment)
Storage:	2–30°C
Tests per kit:	20

Results



Ordering Information

Prod.#	Product description
9719	Reveal for <i>Salmonella Enteritidis</i> (SE) Test System - contains 20 Reveal for SE test devices and 20 disposable transfer pipettes.

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com