

# Neogen's NEW Assay Format

Neogen now offers Forensic Drug Detection Kits in a new and improved format. The new format includes:

- Ready-to-use conjugates
- Calibrators for whole blood, urine and oral fluid
- Extended test kit stability
- Reduced sample dilutions
- Acid stop solution

## Kits currently available in the NEW format:

- Acetaminophen
- Amphetamine Specific<sup>1</sup>
- Amphetamine Ultra<sup>1</sup>
- Barbiturate Group<sup>1</sup>
- Benzodiazepine Group<sup>1</sup>  
(Oxazepam/Clonazepam)
- Buprenorphine
- Carisoprodol
- Cocaine/BZE<sup>1</sup>
- Dextromethorphan
- Fentanyl
- Hydromorphone
- Ketamine
- LSD
- Meperidine
- Methadone
- Methamphetamine/MDMA<sup>1</sup>
- Opiate Group<sup>1</sup>
- Oxycodone/Oxymorphone<sup>1</sup>
- Phencyclidine (PCP)<sup>1</sup>
- Phenytoin
- Propoxyphene
- Sertraline
- THC<sup>1</sup>
- Tramadol
- Tricyclics Group
- Zolpidem

<sup>1</sup>Matrix Specific Calibrators included



944 Nandino Blvd • Lexington, KY 40511  
800/477-8201 USA/Canada or 859/254-1221  
Fax: 859/255-5532 • e-mail: [inform@neogen.com](mailto:inform@neogen.com) • [www.neogen.com](http://www.neogen.com)