

# Colitag

## Presence/Absence (P/A) Water Test Kit




### Intended Use

The Colitag P/A Water Test Kit uses an EPA-approved selective and differential medium to detect total coliforms and *E. coli* in water samples in 16–48 hours.

### Product Specifications

Sensitivity:	1 CFU in a 100 mL sample
Testing time:	16–48 hours
Storage:	2–30°C
Tests per kit:	20, 100 or 250
Approvals:	EPA: Federal Register/vol. 74, No. 216/Tuesday, November 10, 2009/Rules and Regulations

### Results



### Ordering Information

Prod.#	Product description
9850	Colitag P/A water test kits, 20 pack
9851	Colitag P/A water test kits, 100 pack
9852	Colitag P/A water test kits, 250 pack
9855	Sample containers, flexible plastic with flip-top with sodium thiosulfate, 100 pack
9857	Colitag comparator - 100 mL
9859	Sample containers, solid plastic with screw-top with sodium thiosulfate, 100 pack

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)