




Getting Started with Raptor®

Ordering Information

Part#	Description
9680	Raptor Platform
9696	Raptor Solo Platform



Raptor Accessories

Part#	Description
9687	Raptor printer
9689	Raptor keyboard

Extraction Equipment

Part#	Description
9390	Rock-It shaker
9453	Bunn grinder
9401	AgriGrind
9427	Scale 400 g ± 0.1 g
9426	Timer - 3 channel
9475	Reveal® sample cup rack
9330	Mini centrifuge
9447	Graduated Cylinder, plastic

Extraction Consumables

Part#	Description
9428	Sample collection cups
9421	Sample collection tubes
9420	Filter syringes
9172	2.0 mL centrifuge tubes (box 500)
9681	Raptor cartridges (box 50)

Pipettors and Tips

Part#	Description
9860	Pipettor, 100 µL fixed for Q+
9488	Pipettor, 200 µL fixed
9407	Pipette tip rack, 1–200 µL (96 tips)
9410	Pipette tips, 1–200 µL (bag 1000)
9692	Pipettor, 250 µL fixed
9693	Pipettor, 400 µL fixed
9336	Pipettor, 500 µL fixed
9337	Pipettor, 1000 µL fixed
9464	Pipette tips, 100 µL to 1 mL (bag 1000)
9487	Pipette tip rack, 100–1000 µL, sterile (96 tips)
9463	Pipettor, 100–1000 µL adjustable

Pipettors and pipette tips required for Reveal Q+ and Reveal Q+ MAX test kits

Part #	Reveal® Kit Description	100 µL pipette	200 µL pipette	250 µL pipette	400 µL pipette	500 µL pipette	1000 µL pipette	100–1000 µL pipette
		9860	9488	9692	9693	9336	9337	9463
		Pipettors require 1–200 µL pipette tips (9407 & 9410)			Pipettors require 100–1000 µL pipette tips (9464 & 9487)			
8085	Q+ for Aflatoxin	•			•	•		
8385	Q+ for DON	•			•	•	•	
8885	Q+ for Fumonisin		•		•			
8088	Q+ MAX for Aflatoxin			•	•			
8388	Q+ MAX for DON	•			•	•	•	
8888	Q+ MAX for Fumonisin							
8288	Q+ MAX for T-2/HT-2	•			•	•	•	
8188	Q+ MAX for Zearalenone		•	•	•			
8688	Q+ MAX for Ochratoxin	•	•		•			
8488	Q+ MAX for Ergot Alkaloids				•			•

