



Prima[®]-Flex

Product Numbers: See product tables

Presentation

Product Number	Product Description
382368	Prima-Flex Premium Injector 20mL
379566	Prima-Flex Premium Injector 15mL
382351	Prima-Flex Premium Injector 10mL

Product Description

The Prima-Flex is a durable, impact-resistant syringe that quickly and easily converts to a line feed or bottle mount medication application system. The self-contained, one-handed unit is ideal for full mobility and automatic refilling of syringes without the medication returning to the bottle. The Prima-Flex family of applicators features three core barrel sizes: 10mL, 15mL, and 20mL. Each applicator comes standard with metal needle nut, 20mm bottle collar, 30mm bottle collar, 33mm bottle collar, 20mm fast fit, 30mm fast fit, two hose armor springs, soft medical grade tubing with draw off spike, bottle mount configuration with spike, and soft medical grade tubing for bottle mount configuration and is packaged in a premium cardboard box with lubricant packet.



Line Feed Adaptability



Bottle Mount Adaptability

Replacement Parts

Product Number	Product Description
382528	20mL Prima-Flex Replacement Barrel and O-ring
381361	15mL Prima-Flex Replacement Barrel and O-ring
382542	10mL Prima-Flex Replacement Barrel and O-ring
379955	Prima-Flex Replacement Conversion Kit

To learn more, call 800.621.8829 (USA/Canada) or 859.254.1221, or visit neogen.com.

