



PrimaBond™

Product Numbers: See product table below

Key Features

- Use with wooden or rubber blocks, or plastic slip-ons
- Color change confirms proper mixing

Product Description

PrimaBond™ is a fast setting (20-30 seconds) urethane adhesive designed for use in applying wooden or rubber blocks as an adjunct in the treatment of bovine lameness. Product turns green with proper application. PrimaBond is available in a 200 mL tube. Sold in eaches; 12 per case.

Cold Weather Suggestions

- Keep container warm and the material will be easier to work and set at normal times
- Heat with heat gun for 15 seconds
- Keep hoof off ground for an extra 30 seconds

PrimaBond	Product Number
Fast-Set Adhesive 200 mL	SB2012
Mixing Tips, 12/bag	SBT012
Mixing Tips, 100/bag	SBT100
Adhesive Applicator Gun	TA4158

Use with

Wooden Block 1 ct	HC400
Rubber Block 1 ct	RB100

To learn more, call 800.621.8829 (USA/Canada) or 859.254.1221, or visit neogen.com.

