



Ideal[®] Blood Collection Tubes

Product Numbers: BCT-10P, BCT-10P, BCT-5P, BCT-7P, BCT-10EDTA, BCT-4EDTA, BCT-6EDTA, BCT-6SST, BCT-9SST



Description

Our disposable blood collection vacuum tubes are made of shatter-resistant plastic and are latex free.

The tubes are labeled for easy marking, and the interiors are coated with water soluble silicone. Three varieties available. Sterile.



Features

Standard Blood Collection Tubes: Red top tubes contain no additives and are used primarily for serum lab tests.

EDTA Blood Collection Tubes: Lavender top tubes contain the additive EDTA/K2 and are used primarily for CBC lab tests.

SST Blood Collection Tubes: Yellow/gold top tubes contain the additive Polymer Gel. After centrifugation, the gel forms a barrier between the clot and serum.



Standard Blood Collection Tubes

Product Number	Product Description
BCT-10P	Blood Collection Tubes Plain 10 mL Red 100/tray
BCT-3P	Blood Collection Tubes Plain 3 mL Red 100/tray
BCT-5P	Blood Collection Tubes Plain 5 mL Red 100/tray
BCT-7P	Blood Collection Tubes Plain 7 mL Red 100/tray

EDTA Blood Collection Tubes

Product Number	Product Description
BCT-10EDTA	Blood Collection Tubes EDTA 10 mL Pur 100/tray
BCT-4EDTA	Blood Collection Tubes EDTA 4 mL Pur 100/tray
BCT-6EDTA	Blood Collection Tubes EDTA 6 mL Pur 100/tray

SST Blood Collection Tubes

Product Number	Product Description
BCT-6SST	Blood Collection Tubes SST 6 mL Yel 100/tray
BCT-9SST	Blood Collection Tubes SST 9 mL Yel 100/tray

To learn more, call 800.621.8829 (USA/Canada) or 859.254.1221, or visit NEOGEN.com.

