



# Calf Eze™ Calf Puller

Product Number: 3022



## Product Description

Calf Eze™ is engineered for the user that prefers a dual action pull, or zig-zag motion, instead of a single pull instrument.

Calf Eze is superior to ratchet-style pullers that rely on grooved shafts with sharp edges and hard-to-clean grooves. Calf Eze utilizes the same non-slip cam traction system as the Dr. Frank's Calf Puller, but has two pulling points that alternately advance to shift the pull from one leg to the other.

The dual-action jack body is made of special strength aluminum and weighs only 17 lbs. The jack release mechanism allows the user to easily position the jack on the rods and lock it in place. Steel cams prevent slipping and provide longer jack life. A third chain hook is also provided to give a single-pull option.

The Calf Eze is suggested for use with two 45" OB chains (not included; available as item 3102).

## Key Features

- Lightweight - only 17 lbs
- Easy pin release
- Zig-zag pulling points

| Calf Eze Calf Puller        | Product Number |
|-----------------------------|----------------|
| Complete unit               | 3022           |
| Individual Items            | Product Number |
| Jack Assembly with handle   | 3024           |
| 72" Rods (pair)             | 3006           |
| 23" Breechen                | 3008           |
| Web Strap                   | 3009           |
| Calf Eze Cams & Springs Kit | 9785           |
| OB Chain 45"                | 3102           |

To learn more, call 800.621.8829 (USA/Canada) or 859.254.1221, or visit [neogen.com](http://neogen.com).

