

FOR LABORATORY USE

STORE AT 2-8°C

soleris®

Burkholderia Cepacia Complex Supplement

Part Number: S2-BCC-S

20 mL per vial

USAGE INSTRUCTIONS

1. Swirl supplement prior to use.
2. Add 0.2 mL of the S2-BCC-S supplement to the desired PD-109 *Pseudomonas* test vial.
3. Place the remaining S2-BCC-S supplement in the refrigerator.

For more information contact your local NEOGEN representative at 800.234.5333 or visit our website at NEOGEN.com.

