



Raptor[®] Integrated Analysis Platform



Raptor Integrated Analysis Platform

The Raptor® Integrated Analysis Platform is a lateral flow test strip reader with built-in incubation. The system provides an easy way to objectively analyze and store results of Neogen's lateral flow tests. The system controls the timing, temperature and reading of the test strips, allowing users to simply add sample and walk away. The Raptor system protects the integrity of your data by processing and analyzing results without additional operator input.

- Barcoded test strips identify the test type, lot and expiry
- Add sample and walk away
- System controls assay parameters for data integrity
- Connectivity options for easy testing data integration
 - Wi-Fi or USB transfer of results
- Utilizes Neogen Data Manager Software system for easy data analysis

NEOGEN offers two options to meet your needs



Scalability

Three ports for independent processing of three samples



Portability

Battery operated for complete portability

Raptor's patented cartridge

- Provides flexibility for your testing
- Ensures correct volume of sample is tested
- Made with recyclable plastic components



Ordering Information

9680	Raptor Integrated Analysis Platform	8088	Reveal Q+ MAX for Aflatoxin
9696	Raptor Solo Integrated Analysis Platform	8388	Reveal Q+ MAX for DON
9681	Raptor Cartridges	8688	Reveal Q+ MAX for Ochratoxin
9687	Optional Printer	8188	Reveal Q+ MAX for Zearalenone
9689	Optional Keyboard	8288	Reveal Q+ MAX for T-2/HT-2
8085	Reveal Q+ for Aflatoxin	8888	Reveal Q+ MAX for Fumonisin
8385	Reveal Q+ for DON	8488	Reveal Q+ MAX for Ergot Alkaloids
8885	Reveal Q+ for Fumonisin		

Traditional Testing

Other methods of testing require multiple hands-on steps that are potential source of operator error and require close monitoring of the process. All of these touch points have the potential to jeopardize the integrity of your data.

Raptor Platform

The Raptor platform improves the efficiency of your laboratory by eliminating operator error and increasing productivity — just add sample and walk away. The Raptor platform protects the integrity of your data.



