



Aluminum Hub Needles

Product Numbers: See table below



Key Features:

- Aluminum hub for large animal applications requiring a durable needle
- Tri-beveled, anti-coring tips are computer-ground for consistent sharpness
- Anti-friction coating for minimal penetration resistance and animal discomfort
- Stainless steel cannula
- Sterile, non-pyrogenic, auto-clavable components
- Tamper evident cartridges are safe, convenient, easily disposed of and transportable
- Lot number is printed on each cartridge for easy tracking
- ISO color-coded cartridge for easy gauge identification

Needles are available in boxes of 100; master cases contain 10 boxes. Ask about additional size and length availability.

Product Number	Size	ISO Gauge Color
9400AL	14 G x 1" • 2,0 x 25 mm	Green
9401AL	14 G x 1½" • 2,0 x 38 mm	Green
9402AL	14 G x 2" • 2,0 x 51 mm	Green
9403AL	16 G x ¾" • 1,6 x 16 mm	White
9404AL	16 G x ¾" • 1,6 x 19 mm	White
9405AL	16 G x 1" • 1,6 x 25 mm	White
9406AL	16 G x 1½" • 1,6 x 38 mm	White
9407AL	18 G x ⅝" • 1,3 x 16 mm	Pink
9408AL	18 G x ¾" • 1,3 x 19 mm	Pink
9409AL	18 G x 1" • 1,3 x 25 mm	Pink
94010AL	18 G x 1½" • 1,3 x 38 mm	Pink
06-2043B	20 G x ¾" • 0,9 x 19 mm	Yellow
9413AL	20 G x 1" • 0,9 x 25 mm	Yellow
9414AL	20 G x 1½" • 0,9 x 38 mm	Yellow

To learn more, call 800.621.8829 (USA/Canada) or 859.254.1221, or visit NEOGEN.com.

