



Syringe



Prima® Pour-On Drenchers

Product Numbers: See product table below

Product Description

Prima® Pour-On Syringes

These Pour-On applicators (or syringes) allow for adjustable dosage increments with an accurate and easy to use dial-a-dose system. The syringes have a shower head nozzle which avoids wastage from splashing as well as an angled nozzle for easy administration. A T-bar applicator, to more broadly spread liquid down the backline, is also available.

Prima Drenchers

Prima oral drenchers are built at next level standards, suitable for many drenching applications where accuracy and speed are essential. The drenchers have extremely rapid fill time and the dial-a-dose selection allows changes from large to small doses without any product wastage. Accessories include angled, hooked, and floating hook nozzles.

Premium Line Pour-On and Premium Line Drenchers

Packaged with soft medial grade tubing, two hose armor springs, and lube packet. Packaged in a clear heat-sealed bag.

Item Number	Product Description	Size
Pour-On Applicators		
352910	with T-Bar and 8 mm hose barb	30 mL
330901	with Standard nozzle and 8 mm hose barb	30 mL
351098	with Standard nozzle and 12 mm hose barb	30 mL
201058	with Standard nozzle and 12 mm hose barb	50 mL
332585	with Standard nozzle and 12 mm hose barb	70 mL
345745	with Standard nozzle and 12 mm hose barb (various colors)	70 mL
346797	with Standard nozzle and 12 mm hose barb and header card	70 mL
351630	with Standard nozzle and 12 mm hose barb	80 mL
Drenchers		
367037	with 3" metal nozzle	30 mL
366597	with plastic T-Bar nozzle	30 mL
346926	with 8.6" drench hook	30 mL
345103	with 6.6" angled drench hook	70 mL
346568	with 23" remote floating hook	70 mL
365934	with 23" remote floating hook	80 mL
365743	with long metal angled drenched hook	70 mL
365903	with Standard nozzle with 12 mm hose barb	80 mL



Accessories

Prima[®] Pour-On/Drenchers Accessories

Product Numbers: See product table below

Product Description

Additional accessories are available for all of the Prima Pour-Ons and Drenchers providing users with versatility to meet their needs. Quality replacement parts help keep your Prima syringes serviceable for years to come, reducing your operational costs.

Product Image	Item Number	Product Description	Size	Angle
Metal Nozzles				
	342737	Angled Hooked Metal Nozzle	8.6" x 5.9" (21.8 cm x 14.9 cm)	45°
	346957	Angled Hooked Metal Nozzle	6.6" x 4.3" (16.8 cm x 10.9 cm)	35°
	378064	Metal Nozzle	3"	n/a
	375285	Metal Nozzle	11"	35°
	345110	Floating Hook - Remote Drench Hook	23" x 8" (58.4 cm x 20.3 cm)	45°
	351494	Plastic Connector for Floating Hook	n/a	n/a
Caps				
	341044	Vented Cap	43 mm	
	372970	Vented Cap and Dip Tube	5 L/20 L (28 mm/63 mm)	
	372963	Inverted Cap	5 L (28 mm)	
	349066	Vented Cap and Dip Tube	38 mm	
Replacement Barrel Kits				
	353627	Replacement Barrel with O-Ring, Felt, Piston	30 ml	
	346179	Replacement Barrel with O-Ring, Felt, Piston	70 ml	
Plastic Nozzles				
	352507	T-Bar Nozzle		
	365927	Fan Tip Nozzle		

To learn more, call 800.621.8829 (USA/Canada) or 859.254.1221, or visit NEOGEN.com.

