



Repeater Syringes



Prima-Shot™

Product Numbers: See product table below

Key Features

- Amber-colored barrel
- 10 fixed-dose settings
- Durable rugged construction for repetitive use
- Ergonomic design is easy on hands

Product Description

The Prima-Shot™ 50 mL repeater syringe adjusts from 0.5 mL to 5 mL in 0.5 mL increments, and features an easy to grip dose adjuster dial to make dosing changes easy. Our 25 mL syringe adjusts from 0.25 to 2.5 mL in 0.25 increments. Features a heavy duty Tritan™ Barrel that is BPA-free, highly chemical and scratch-resistant, and has chemical compatibility that surpasses both polycarbonate and polypropylene. Includes colored pull knob inserts that can quickly be exchanged to easily distinguish between medications.

Purpose

NEOGEN's Prima-Shot is engineered to be the most accurate repeater syringe in the animal health industry. The patented design presents a unique rotating pull knob to conveniently fill and engage the syringe. Designed with the end user in mind, the ergonomic handles are sure to make any job run smooth.

To learn more, call 800.621.8829 (USA/Canada) or 859.254.1221, or visit NEOGEN.com.

Prima-Shot	Product Number
Repeater 50 mL w/ amber barrel	367723
Repeater 25 mL w/ amber barrel	377210

Replacement Parts

Metal needle nut	330840
Replacement pull knob inserts (4/bag)	374158
Prima-Shot 50 mL amber replacement barrel w/ o-ring	374172
Replacement Prima-Shot 50 mL o-rings (12/bg)	374196
Prima-Shot 25 mL amber replacement barrel w/ o-ring	377883
Replacement Prima-Shot 25 mL o-rings (12/bg)	377890





Repeater Syringes



Prima-Shot™

Product Numbers: See product table below

Key Features

- Clear barrel
- 10 fixed-dose settings
- Durable rugged construction for repetitive use
- Ergonomic design is easy on hands

Product Description

The Prima-Shot™ 50 mL repeater syringe adjusts from 0.5 mL to 5 mL in 0.5 mL increments, and features an easy to grip dose adjuster dial to make dosing changes easy. Our 25 mL syringe adjusts from 0.25 to 2.5 mL in 0.25 increments. Features a heavy duty Tritan™ Barrel that is BPA-free, highly chemical and scratch-resistant, and has chemical compatibility that surpasses both polycarbonate and polypropylene. Includes colored pull knob inserts that can quickly be exchanged to easily distinguish between medications.

Purpose

NEOGEN's Prima-Shot is engineered to be the most accurate repeater syringe in the animal health industry. The patented design presents a unique rotating pull knob to conveniently fill and engage the syringe. Designed with the end user in mind, the ergonomic handles are sure to make any job run smooth.

To learn more, call 800.621.8829 (USA/Canada) or 859.254.1221, or visit NEOGEN.com.

Prima-Shot	Product Number
Repeater 50 mL w/ clear barrel	376428
Repeater 25 mL w/ clear barrel	378002

Replacement Parts

Metal needle nut	330840
Replacement pull knob inserts (4/bag)	374158
Prima-Shot 50 mL clear replacement barrel w/ o-ring	377975
Replacement Prima-Shot 50 mL o-rings (12/bg)	374196
Prima-Shot 25 mL clear replacement barrel w/ o-ring	377982
Replacement Prima-Shot 25 mL o-rings (12/bg)	377890

