


Schedule of Accreditation

issued by

United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

 <p>UKAS TESTING</p> <p>8829</p> <p>Accredited to ISO/IEC 17025:2017</p>	<p>Neogen Europe Limited</p> <p>Issue No: 010 Issue date: 08 February 2021</p>	
	<p>1 Quest Park Moss Hall Road Heywood Lancashire BL9 7JJ</p>	<p>Contact: Samantha Cowley Tel: +44 (0)161 820 3833 Fax: +44 (0)161 820 5383 E-Mail: SCowley@neogen.com</p>
<p>Testing performed at the above address only</p>		

DETAIL OF ACCREDITATION

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>MICROBIOLOGICAL CULTURE MEDIA Ready to use solid and liquid media, including diluents</p>	<p>Performance testing of physical and microbiological parameters to support production quality of finished items for sale</p> <p><u>Microbial Performance Tests</u></p>	<p>Documented In-house methods in accordance with BS EN ISO 11133:2014+A2:2020 and in conjunction with Product Specifications</p>
<p>Solid media</p>	<p>Evaluation of growth promotion on solid media (productivity and selectivity)</p>	<p>1) LM6.9 Qualitative performance testing for solid medium</p> <p>2) LM5.13 Quantitative performance testing for solid medium</p>
<p>Liquid media</p>	<p>Evaluation of growth promotion in liquid media (productivity and selectivity)</p>	<p>LM4.3 Qualitative performance testing for liquid media</p>
<p>Diluents</p>	<p>Maintenance of viability</p> <p><u>Physical Tests</u></p>	<p>LM4.5 Quantitative evaluation of viability maintenance in diluents over a holding time of 45 mins</p>
<p>Solid media</p>	<p>pH</p> <p>Fill volume weight check</p>	<p>LM3.5 pH using a pH meter</p> <p>LM3.4 Fill volume weight check for solid medium</p>



8829

Accredited to
ISO/IEC 17025:2005

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
21 - 47 High Street, Feltham, Middlesex, TW13 4UN, UK

Neogen Europe Limited

Issue No: 010 Issue date: 08 February 2021

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>MICROBIOLOGICAL CULTURE MEDIA Ready to use solid and liquid media, including diluents</p> <p>Solid media</p> <p>Liquid media</p>	<p>Performance testing of physical and microbiological parameters to support production quality of finished items for sale</p> <p><u>Physical Tests (Cont'd)</u></p> <p>Contamination check at 25°C and 37°C</p> <p>pH</p> <p>Contamination check at 25°C and 37°C</p>	<p>Documented In-house methods in accordance with BS EN ISO 11133:2014+A2:2020 and in conjunction with Product Specifications</p> <p>LM3.6 Sterility test method for solid medium</p> <p>LM4.2 pH using a pH meter</p> <p>LM4.1 Assessment of growth on TSA following incubation</p>
END		